



# Quantitative Real-time Polymerase Chain Reaction (Coupled with High-resolution Melting) for Simultaneous Detection and Quantification of Four *BCR::ABL1* Isoforms

T. Zhou<sup>1</sup>, A. Mohamed<sup>2</sup>, B. Legendre<sup>2</sup>, H. Newell<sup>2</sup>, K.S.J. Elenitoba-Johnson<sup>1</sup>

1. Department of Pathology and Laboratory Medicine, Memorial Sloan Kettering Cancer Center, New York, NY, USA

2. Precipio Inc, New Haven, Connecticut, USA.



## INTRODUCTION

The diagnosis and management of chronic myeloid leukemia (CML) critically rely on detecting and serially quantifying *BCR::ABL1*.

Currently, no clinical assay simultaneously detects multiple *BCR::ABL1* isoforms; each isoform requires a separate analysis, which is labor-intensive and cost-inefficient.

This technical constraint results in several important limitations in clinical practice:

- Baseline limits: testing is typically limited to p210/ p190, which may miss p230 or atypical isoforms as primary drivers when co-expressed with low-level p210/p190.
- Monitoring gap: standard monitoring tracks only baseline transcripts, potentially overlooking emerging novel isoforms during evolution, falsely suggesting disease stability and delaying therapeutic adjustments.
- Coverage gap: the prevalence and significance of e13a3/e14a3 (p203)—anecdotally linked to asciminib resistance—remain unclear.

## AIM

- To develop a single multiplex assay that simultaneously detects and quantifies four clinically relevant *BCR::ABL1* transcripts—e1a2/e1a3 (p190), e13a2/e14a2 (p210), e19a2/e19a3 (p230), and e13a3/e14a3 (p203).
- To evaluate its analytical performance and clinical utility.

## METHOD

We developed the **BloodHound** assay, employing multiplex RT-qPCR. This assay achieves a limit of detection of 0.001%, 99.9% sensitivity, and > 99% specificity, when validated against droplet digital PCR and Sanger sequencing.

Using **BloodHound**, we characterized the distribution and co-expression of the four transcripts in 895 samples (**Table 1**), from patients with suspected (n = 566), established (n = 296), relapsed (n = 32) CML, or unknown status (n = 1).

Sanger sequencing and Qiagen ipsogen IS assay were used for orthogonal confirmation.

## RESULTS

**Overall *BCR::ABL1* detection:** 187/895 samples (20.9%)

- Suspected CML: 3.5% (20/566)
- Established CML under monitoring: 33.1% (98/296)
- Relapsed disease: 100% (32/32)


### Isoform distribution

- p210 alone: 161/187 (86.1%)
- p190 alone: 2/187 (1.1%)
- p230 alone: 1/187 (0.5%)
- P203: not detected.
- p190 and p210 co-expression: 23/187, 12.3%) was more frequent than historically appreciated, particularly in baseline diagnostic specimens (25.0%).
- When co-expressed, p210 levels were markedly higher than p190 (median p210:p190 ratio = 2,152).

### Orthogonal confirmation:

- Qualitative results showed 100% concordance with Sanger sequencing.
- Quantitative values correlated tightly with the standardized Qiagen ipsogen IS assay (r = 0.998; median absolute difference = 0.02%, **Fig. 1**).

**Figure 1.** Correlation of p210 transcript levels between the BloodHound™ assay and the parallel Qiagen assay.



**Table 1. Specimen types and disease status**

Characteristics	Values
Total specimens tested; N	895
Specimen types; N (%)	
BM	49 (5.5%)
PB	846 (94.5%)
Overall disease status; N (%)	
Suspected/baseline CML	566 (63.2%)
Post-therapy	296 (33.1%)
Relapse	32 (3.6%)
Unknown	1 (0.1%)

**Table 2. Distribution of *BCR::ABL1* transcripts**

<i>BCR::ABL1</i> positivity	n/N (%)
Total positive specimens	187/895 (20.9%)
p190 alone	2/187 (1.1%)
p210 alone	161/187 (86.1%)
p230 alone	1/187 (0.5%)
p203 alone	0
p190 + p210 co-expression	23/187 (12.3%)
All other co-expression	0
Suspected/baseline CML	20/566 (3.5%)
p190 alone	0
p210 alone	14/20 (0.0%)
p230 alone	1/20 (5.0%)
p190 + p210 co-expression	5/20 (25.0%)
Post-therapy	135/296 (45.6%)
p190 alone	2/135 (1.5%)
p210 alone	118/135 (87.4%)
p230 alone	0
p190 + p210 co-expression	15/135 (11.1%)
Relapse	32/32 (100.0%)
p190 alone	0
p210 alone	29/32 (90.6%)
p230 alone	0
p190 + p210 co-expression	3/32 (9.4%)

## CONCLUSIONS

- BloodHound simultaneously detects and quantifies four *BCR::ABL1* isoforms in a single reaction, addressing major gaps in CML testing.
- The assay shows robust clinical performance and high analytical concordance with standard methods and across diagnostic, monitoring, and relapse settings.
- Co-expression of p190 and p210 is more common than previously reported
- BloodHound has the potential to significantly streamline clinical workflows, enhance diagnostic precision, and provide deeper insights into CML clonal dynamics.

## REFERENCES

- Rowley JD. Letter: A new consistent chromosomal abnormality in chronic myelogenous leukemia identified by quinacrine fluorescence and Giemsa staining. *Nature*. 1973;243(5405):290-3
- Zhou T et al. Chronic Myeloid Leukemia: Beyond BCR-ABL1. *Curr Hematol Malig Rep*. 2018;13(6):435-45
- Baccarani M et al. The proportion of different BCR-ABL1 transcript types in chronic myeloid leukemia. An international overview. *Leukemia*. 2019;33(5):1173-83
- Leske IB et al. The e13a3 (b2a3) and e14a3 (b3a3) *BCR::ABL1* isoforms are resistant to asciminib. *Leukemia* 2024;38(9):2041-5

## Technical and Clinical Advantages of Bloodhound

- Single-reaction multiplexing
- Pre-plated panel with embedded primers and controls, simplifying implementation and reducing handling errors
- High analytical performance
- Compatible with standard laboratory equipment
- Robust intra-laboratory, inter-laboratory, and inter-primer set standardization
- Potentially interlaboratory standardization of non-p210 quantification

## ACKNOWLEDGEMENTS

The authors gratefully acknowledge the support and collaboration of Dr. Ting and Dr. Elenitoba-Johnson of the Memorial-Sloan Kettering Cancer Center, for their substantive scientific, clinical and editorial insights that made this publication possible.

## CONTACT INFORMATION



Precipio, Inc  
WWW.Precipiodx.com  
[Contact@precipiodx.com](mailto:Contact@precipiodx.com)  
203-787-7888

