

Precipio's Hematologic Malignancies Panel focuses on clinically actionable mutations of interest. Unlike other NGS panels with varied limits of detection across mutation targets, Precipio's panel provides deep, uniform coverage:

1% limit of detection for all genes




Designed to support our Omnia methodology to deliver the most accurate diagnosis for your patients.

Panel Genes

| | | | | |
|---------|--------|---------|--------|---------|
| ABL1 | BRAF | CEBPA• | ETV6• | HRAS |
| KDM6A• | NPM1 | PTEN | SMC1A | TP53• |
| ANKRD26 | CALR | CSF3R | EZH2• | IDH1 |
| KIT | NRAS | PTPN11 | SMC3 | U2AF1 |
| ASXL1 | CBL | CUX1 | FLT3 | IDH2 |
| KMT2A | PDGFRA | RAD21• | SRSF2 | WT1 |
| ATRX | CBLB | DDX41 | GATA1 | IKZF1• |
| KRAS | PHF6• | RUNX1• | STAG1 | ZRSR2• |
| BCOR• | CBLC | DNMT3A• | GATA2• | JAK2 |
| MPL | PIGA• | SETBP1 | STAG• | BCORL1• |
| CDKN2A• | ETNK1 | GNAS | JAK3 | NF1 |
| PPMID | SF3B1 | TET2• | | |

• Full CDS coverage in gene indicated by

Specimen Requirements:

| Specimen Type | Requirements | Transportation Conditions |
|--|--|---|
|  Bone Marrow | 2 mL EDTA Tube preferred | <ul style="list-style-type: none"> Cold pack <ul style="list-style-type: none"> No direct contact with specimen. Ship as soon as sample is collected. Overnight delivery |
|  Peripheral Blood | 5 mL EDTA Tube preferred | |
|  Formalin Fixed Paraffin Embedded Tissue (FFPE) | 1 H&E slide and 6-10 unstained slides, 5-7 microns of BM clot or tissue fixed with 10% NBF fixative. Or, the BM clot FFPE block can be sent for sectioning in our lab. | <ul style="list-style-type: none"> Block packed in bubble-wrap plastic Slides placed in plastic slide holder |

Documentation:

Completed & signed Precipio requisition
 Patient clinical history
 Prior pathology reports
 US customs form (Outside US)
Climate control not required

Ship by registered courier to:

Precipio Inc.
 Attn: Accessioning Department
 4 Science Park | New Haven, CT 06511, USA

For more information please contact:

contact@precipiodx.com | www.precipiodx.com | Tel: +1.203.787.7888

rev 2307