

NGS Hematologic Malignancies Panel

Targeted panel with uniform, deep sequencing coverage

Precipio's Hematologic Malignancies Panel focuses on clinically actionable mutations of interest. Unlike other NGS panels with varied limits of detection across mutation targets, Precipio's panel provides deep, uniform coverage:

TAT 5 Business Days



**Prior Authorization
Support**





1% limit of detection for all genes

Panel Genes

ABL1	ANKRD26	ASXL1	ATRX	BCOR•
BCORL1•	BRAF	CALR	CBL	CBLB
CBLC	CDKN2A•	CEBPA•	CSF3R	CUX1
DDX41	DNMT3A•	ETNK1	ETV6•	EZH2•
FLT3	GATA1	GATA2•	GNAS	HRAS
IDH1	IDH2	IKZF1•	JAK2	JAK3
KDM6A•	KIT	KMT2A	KRAS	MPL
NF1	NPM1	NRAS	PDGFRA	PHF6•
PIGA•	PPMID	PTEN	PTPN11	RAD21•
RUNX1•	SETBP1	SF3B1	SMC1A	SMC3
SRSF2	STAG•	STAG1	TET2•	TP53•
U2AF1	WT1	ZRSR2•		

• Indicates full CDS coverage

Specimen Requirements:

Specimen Type	Requirements	Transportation Conditions
 Bone Marrow Aspirate	2 mL EDTA Tube Only	<ul style="list-style-type: none">• Cold pack<ul style="list-style-type: none">◦ No direct contact with specimen.• Ship as soon as sample is collected.• Overnight delivery
 Peripheral Blood	5 mL EDTA Tube Only	

Documentation:

Completed & signed Precipio requisition
Patient clinical history
Prior pathology reports
US customs form (Outside US)
Climate control not required

Ship by registered courier to:

Precipio Inc.
Attn: Accessioning Department
4 Science Park | New Haven, CT 06511, USA